

# RNA-Seq Technical Specifications




## RNA Sequencing Services

- Standard RNA-Seq
- Strand-Specific RNA-Seq
- Small RNA-Seq
- Ultra-Low Input RNA-Seq
- Single-Cell RNA-Seq\*
- Long-Read RNA-Seq using PacBio or Oxford Nanopore Technologies\*
- Spatial Transcriptomics\*
- CLIA RNA-Seq\*

\*See [genewiz.com](https://www.genewiz.com) for technical specifications.

## RNA Sequencing Workflow



### 1 Experimental Design

We offer resources to help you find the best NGS solution and experimental design for your project.



Interactive NGS Solution Selection Tool:  
[genewiz.com/ngs](https://www.genewiz.com/ngs)



Contact us for a **free technical consultation** with a Ph.D.-level scientist  
[web.genewiz.com/ngs-inquiry](https://www.web.genewiz.com/ngs-inquiry)

\*Other sample types accepted. View [Sample Submission Guidelines](#) for details.

†Please inquire about submitting lower inputs.


\*Contact us about our RNA Stabilization Tubes to ship RNA samples at ambient temperature.

### 2 Sample Preparation

Our extraction solutions are compatible with over 30 standard and hundreds of custom sample types to accommodate your starting material. The below table is by no means an exhaustive list.

| Sample Type*                   | Minimum Amount†   | Recommended Amount   |
|--------------------------------|---|--|
| <b>Total RNA<sup>‡</sup></b>   | Standard: 150 ng, 25 µL, ≥6 ng/µL<br>Ultra-Low: ≥50 pg, 3-5 µL, ≥10 pg/µL | Standard: ≥500 ng, 25-50 µL, 20-50 ng/µL<br>Ultra-Low: ≥10ng, 3-5µL, ≥2ng/µL |
| <b>Eukaryotic cell pellet</b>  | 1×10 <sup>4</sup> cells   | ≥1×10 <sup>6</sup> cells   |
| <b>Prokaryotic cell pellet</b> | 1×10 <sup>6</sup> cells   | ≥1×10 <sup>8</sup> cells   |
| <b>Frozen tissue</b>           | 2 mg  | 30 mg  |
| <b>FFPE</b>                    | 2 curls, 5µm-20 µm, >150 mm <sup>2</sup>                                  | ≥4 curls, 5µm-20 µm, >150 mm <sup>2</sup>                                    |

# RNA-Seq Technical Specifications



## 3 Library Preparation

| RNA-Seq Service            | Target RNA                      | RNA Selection Method  |
|----------------------------|---------------------------------|---|
| Standard & Strand-Specific | mRNA (eukaryotic)               | Poly(A) selection   |
|                            | mRNA + lncRNA                   | rRNA depletion  |
| Small                      | Small RNA (miRNA, siRNA, piRNA) | Size fractionation with adapter ligation to 5' phosphate      |
| Ultra-Low Input            | mRNA (eukaryotic)               | Poly(A) selection with enrichment for full-length transcripts |

## 5 Data Analysis

| RNA-Seq Service                                | Standard Analysis Package  | Additional Analysis Options   |
|--|--|---|
| Standard<br>Strand-Specific<br>Ultra-Low Input | <ul style="list-style-type: none"> <li>Trimming</li> <li>Mapping</li> <li>Differential gene expression</li> </ul>                              | <ul style="list-style-type: none"> <li>Gene fusion discovery</li> <li>RNA SNP/INDEL detection</li> <li>Novel transcript discovery</li> <li><i>De novo</i> transcriptome assembly</li> </ul> |
| Small  | <ul style="list-style-type: none"> <li>Trimming</li> <li>Mapping</li> <li>Differential gene expression</li> <li>Small RNA discovery</li> </ul> |   |

## 4 Sequencing

|               |  |
|---------------|--|
| Platform      | Illumina® NovaSeq™                       |
| Configuration | 2x150 bp                                 |
| Depth         | Customizable to your project needs*      |
| Data Quality  | Guaranteed ≥85% bases with Q30 or higher |

\*Generally, we recommend 5-10 million read pairs per sample for small genomes (e.g. bacteria) and 20-30 million read pairs per sample for large genomes (e.g. human, mouse). Medium genomes often depend on the project, but 15-20 million read pairs per sample is typically sufficient. For *de novo* transcriptome assembly projects, we recommend 100 million read pairs per sample.

## 6 Project Delivery

| Deliverables for All Projects   | Optional Deliverables  |
|---|--|
| <ul style="list-style-type: none"> <li>Sample quality control report</li> <li>Raw data (FASTQ files)</li> </ul> | <ul style="list-style-type: none"> <li>Aligned data (BAM file)</li> <li>Hit counts (TXT file)</li> <li>DGE results (CSV file)</li> <li>GO enrichment analysis (CSV file)</li> <li>Differential splicing analysis (DEXSeq report)</li> <li>De-multiplexed, aggregated Picard BAM file with summary metrics</li> </ul> |


### DATA DELIVERY OPTIONS



SFTP



Customer Cloud Account



External Hard Drive